

Membrane Tacks

Membrane Tacks are used for the easy and secure attachment of regenerative membranes to the bone during guided bone regeneration procedures.

Features & benefits

Flexible – A solution for soft to hard bone.


Easy handling – Fast and easy membrane fastening and removal. For Tacks use the Tack positioning tool.

Easy storage – Cassette for Tacks.

Facts & figures

Packaging – 10-pack, delivered sterile in a blister. tacks are distributed in three separate compartments in the blister for easy selection of appropriate amount for the surgery. Single use.

Material – Titanium alloy grade 5.



Assortment

Art. No.	Description
64036	Membrane Tack 3 mm – 10 pcs
64031	Tack Positioning Instrument
64034	Cassette

All products might not be available in all markets.

See instructions for use for all procedure details.